Cellu M6° Endermolab S

























INTEGRATED LPG COACH

- → Adjustable color touch screen.
- → Step-by-step preconfigured care protocols to visually follow the progress of Endermolift and Lipomassage sessions.
- \rightarrow 180° screen rotating mobilty.

BODY MODULE

New Ergodrive care head

- $\rightarrow \ \, \mbox{\bf Unprecedented depth of tissue grabbing} \\ \mbox{ for greater effectiveness.}$
- → Independent motorized rollers (LPG patents) for high-precision treatment.
- → Control screen, integrated commands, and rotating head for greater usability.
- → 4 rotations of rollers for personalized results (Slim down, Firm, Resculpt, Smooth).

FACE MODULE

3 new Lift care heads.

- → More dynamic cellular stimulation.
- → Adjustable lift for each skin type.
- → 3 intensity settings for high-precision redensifying (Clarify, Replenish, Firm).



ADAPTABLE CONNECTIVITY

- → Plug & Care technology for greater flexibility.
- $\rightarrow \ \, \text{USB port for } \mathbf{regular \, LPG \, update} \\ \mathbf{downloads} \, (\text{new protocols}).$


























Cellu M6° Endermolab S Body and face care

High-technology slimming and anti-ageing

PREMIUM TREATMENT TABLE

- → Adjustable height and reclining back-rest.
- ightarrow White leather and chrome finish.
- → Prestige comfort.



Endermologie is an exclusive technique of **cellular stimulation** by mechanotransduction. Created by LPG in 1986, it enables 100% natural care of fat storage and ageing skin.



ERGODRIVE HEAD

The driving seat of the system which accommodates 3 different KM Roll.



KM ROLL

Exclusive LPG patent: independently motorized rollers.

KM80: (NEW) specifically for large

areas (fat reduction and smoothing).

KM70: Firming and anti-ageing.

KM50: specifically for narrow areas (fat reduction & smoothing).



LIFT HEADS

TML 10: Lips and eye contours.


TML 20: Entire face.

TML 30: Chest and neckline.



AUXILIARY HEADS

For access to smaller areas.



Tel.: +33 (0)4 92 38 39 00 Fax: +33 (0)4 92 96 09 65



Cellu M6[®] Endermolab 🗓

























INTEGRATED LPG COACH

- → Adjustable color touch screen.
- \rightarrow Pre-configured treatment protocols.
- \rightarrow 180° screen rotating mobilty.

THERAPEUTIC MODULE

New TR50 treatment head

- ightarrow Specially designed for therapeutic applications.
- → Independent motorized rollers (LPG patents)
- $\rightarrow \ 3 \ \text{settings for high-precision } \ \text{tissular grabbing} \\ \ \text{(deep or superficial)}.$
- $\,\,\rightarrow\,\,$ Control screen and integrated commands.

BODY MODULE

New Ergodrive treatment head

- $\begin{array}{l} \rightarrow & \textbf{Unprecedented depth of tissular grabbing} \\ & \text{for greater effectiveness.} \end{array}$
- → Independent motorized rollers (LPG patents) for high-precision treatment.
- → Control screen, integrated commands and rotating head for greater usability.
- → 4 rotations of rollers for personalized results (Slim down, Firm, Resculpt, Smooth).



FACE MODULE

3 new Lift treatment heads

- $\rightarrow \ \ \text{More dynamic cellular stimulation.}$
- \rightarrow Adjustable lift for each skin type.
- → 3 intensities for high-precision redensifying (Clarify, Replenish, Firm).

























High-technology therapeutics and aesthetics

■ PREMIUM TREATMENT TABLE

- → Adjustable height and reclining back-rest.
- $\,\, o\,$ White leather and chrome finish.
- → Prestige comfort.

ADAPTABLE CONNECTIVITY

- → Plug & Care technology for greater flexibility.
- ightarrow USB port for **regular LPG update downloads** (new protocols).







Endermologie is an exclusive technique of **cellular** stimulation by mechanotransduction. Created by LPG in 1986, it enables a **100% natural** treatment of **connective** tissue transformations for beauty or therapeutic purposes.



ERGODRIVE HEAD

The driving seat of the system which accommodates 3 different KM Roll.



KM ROLL

Exclusive LPG patent: independently motorized rollers.

KM80: (NEW) specifically for large areas (fat reduction and smoothing).

KM70: Firming and anti-ageing.

KM50: specifically for narrow areas (fat reduction & smoothing).



TR50

Dedicated to therapeutic applications: fibrosis, oedema, inflammation, pre- and post-liposuction etc.



LIFT HEADS

Therapeutic (oedema, fibrosis) and aesthetic (anti-ageing) face treatment.

TML 10: Lips and eye contours.

TML 20: Entire face.

TML 30: Chest and neckline.



AUXILIARY HEADS

For access to smaller and delicate areas (lesions, scars etc).



MICRO-NOZZLES AND MICRO-HEADS

For therapeutic micro-targeted action (fingers for example).



Tel.: +33 (0)4 92 38 39 00 Fax: +33 (0)4 92 96 09 65 www.lpgsystems.com

